HAYER & BOECKER

DIE DRAHTWEBER

HAYER EML 200 digital plus

**ANALYSENSIEBMASCHINE – TEST SIEVE SHAKER –
TAMISEUSE DE LABORATOIRE**

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1. Introduction

Only install and start the test sieve shaker after having read the following operating instructions.

The operating instructions give important information about the secure, appropriate and economic operation of the test sieve shaker. They should be respected in order to avoid any dangers as well as to maximize the reliability and durability of the machine.

The safety instructions and regulations of prevention of accidents, as well as the different protective recommendations of the technical associations and employer's liability insurance associations, have to be strictly observed when using the machine.

**All safety instructions within these operating instructions have to be followed!
Always keep the operating instructions with the machine.**

If used properly, the construction of the machine guarantees that it runs and can be operated without any risk to persons.

The guarantee does not apply to damages caused by improper use due to a non-respect of the operating instructions.

All pictures and drawings in these operating instructions are general illustrations. They do not reflect the construction in detail.

The technical data and dimensions are without obligation, claims cannot be derived.

We reserve technical changes of the machine described below, without changing the operating instructions.

All rights on these operating instructions, the added drawings and other documents remain to HAVER & BOECKER. So are the rights on disposition, copying and passing on, including trademark registering.

2. Transport

The HAVER-Test Sieve Shaker EML 200 digital plus is delivered completely assembled, except the guide rods.

If not otherwise agreed, it is packed according to the packing directives „HPE“, as defined by „Bundesverband Holzmittel, Paletten, Exportverpackungen e.V.“ (the national association of wooden instruments, pallets, export packaging) and „Verein Deutscher Maschinenbauanstalten“ (the society of German engineering shops)






Observe the current valid national and international directives for accident prevention. Damage caused by improper transport will invalidate any guarantee or replacement.

Warning

Improper lifting may cause bodily harm or material damage. Only lift the machine with the suitable equipment and by qualified persons!

3. Quantity delivered
3.1. EML 200 digital plus


<i>EML 200 digital plus T</i>	<i>EML 200 digital plus N</i>
 <p data-bbox="300 566 746 600"><i>Basic machine EML 200 digital plus T</i></p>	 <p data-bbox="973 566 1420 600"><i>Basic machine EML 200 digital plus N</i></p>
 <p data-bbox="368 804 687 837"><i>cover with inspection glass</i></p>	 <p data-bbox="979 804 1414 869"><i>cover with inspection glass and wide spreading spray diffuser</i></p>
 <p data-bbox="261 1178 799 1243"><i>1 pair of guide rods, nuts with handlebars and counter nuts</i></p>	 <p data-bbox="935 1178 1473 1243"><i>1 pair of guide rods, nuts with handlebars and counter nuts</i></p>
 <p data-bbox="331 1447 715 1512"><i>1 spanner SW 27 DIN 894 WGB to fix the guide rods</i></p>	 <p data-bbox="1005 1447 1388 1512"><i>1 spanner SW 27 DIN 894 WGB to fix the guide rods</i></p>
 <p data-bbox="225 1760 836 1825"><i>1 connection cable, 1,5m long, with EURO-plug and appliance plug</i></p>	 <p data-bbox="879 1760 1528 1848"><i>Separate control gear with special plug type, holding device with plugs and screws for wall fixing and 3 lin.m PVC-water hose</i></p>

Only connect the control gear delivered with the test sieve shaker to the drive!

3.2. Clamping systems (according to model)

3.2.1. Classic

1. Nut with handlebars
2. Sieve cover
3. Guide rod


Elevation adjustment:

By turning the nut it moves up and down the guide rod.

Tensioning:

The sieve cover has to be installed and both nuts must be screwed down **simultaneously**.

3.2.2. HAVER TwinNut

1. Nut HAVER TwinNut
2. Guide rod
3. Sieve cover

Elevation adjustment:

Pull apart the two parts of the nut and slip it on the guide rod up to the sieve cover. Reclose it.

Tensioning:

Tighten the two HAVER TwinNuts **simultaneously**.

4. Technical data for the EML 200 digital plus

Test sieve shaker with electromagnetic drive and three-dimensional sieving action for dry and wet sieving.

Operation voltage: 230 Volt or 110 Volt with transformer

Frequency: 50 Hertz / 60 Hertz

Power consumption: approx. 410 Volt-ampere with full charge

Control: **all-electronic**

Weight: 34 kg, without test sieves

max. weight of full sieve tower: 8,7 kg


Dimensions: 345 x 285 x 930 mm (Dept x Width x Height)

Emitted sound: The emitted sound of the machine in action has been determined with a comparative sound measurement according to DIN 45635-01 K11.

The equivalent continuous sound level A is: $L_{pAeq} \leq 70$ (dB)

The intensity of sound can be reduced when using HAVER Sound absorbing Cupboards.

5. Description EML 200 digital plus T (dry sieving)


1. Nut with handlebars
2. Inspection glass
3. Cover
4. Rubber buffer
5. Test sieve set
6. Guide rod
7. Oscillation magnet
8. Rubber profile for vibration body
9. Vibration body
10. Metal feet rubberized
11. Leaf spring
12. Base plate
13. Front panel with control board
14. Rubber ring
16. Sieve pan
17. Plug connection with main switch and fuse box
18. Counter nut

Never use the EML 200 digital plus T in damp rooms! Use instead the EML 200 digital plus N for wet sieving, with separate control gear!

6. Description EML 200 digital plus N (wet sieving)

1. Nut with handlebars
2. Inspection glass
3. Cover
4. Rubber buffer
5. Test sieve set
6. Guide rod
7. Oscillation magnet
8. Rubber profile for vibration body
9. Vibration body
10. Metal feet rubberized
11. Leaf spring
12. Base plate
13. Separate control unit
14. Rubber ring
15. Wide spreading spray diffuser
16. Sieve pan with outlet
17. Plate with cable with screw fitting
18. Counter nut

With the following parts the EML 200 digital plus N is also suitable for dry sieving:

- Sieve cover with inspection glass, Art. No.: 00561521
- Sieve pan without outlet (see annex, table)

7. PC – Connection

For data transfer the EML 200 can be connected with a personal computer by its serial interface. A RS 232 cable is not included in delivery.

The serial interface of the EML 200 digital plus for dry sieving is located at the backside of the drive. The serial interface of the EML 200 digital plus for wet sieving is located at the separate control gear.

This enables the communication to report the analysis with the Haver software analysis or any other laboratory information and management system (LIMS).

Serial interface at the bottom of the separate control gear.

8. Application references


The machine has been constructed according to the actual technical standards and the recognized safety specifications. Although, there might be dangers for the user or third persons or impairments of the machine and other goods.

Only use the machine in perfect technical constitution. Follow the operating instructions conscious of safety and danger. Repair defects at once by authorized and qualified persons!

The machine has exclusively been developed for grain sieving to determine particle distribution. Any other or further application is not as agreed upon. The manufacturer/supplier will assume no responsibility for damages as a result of improper use. The risk is on the user's side.


9. Assembly

1. Put the machine on a **stable, vibration free table. Do not use a rubber support!**
2. Screw the guide rods (3) into the drilled holes of the vibration plate (1) and drive in the counter nuts (2) by the key.

Minimum of three test sieves (height 50 mm) and one sieve pan is necessary when using guide rods with lengths of 660 mm

EML 200 digital plus T


1. Put on the test sieve set (1) and fill the sieving material in the upper test sieve.
2. Put on the cover (2) and fix it **simultaneously** with the two tension nuts (3).

Please notice item 10, Test sieving according to ISO 2591-1!

Connect the mains plug, pay attention to the power supply!

EML 200 digital plus N

1. Assembly of wide spreading spray diffuser: put the wide spreading spray diffuser (6) through the inspection glass (6); fit the sealing ring (7) and the hose clip (8); **tighten**
2. Put on the sieve pan with outlet (1), the test sieve set (2) and fill the sieving material in the upper test sieve.
3. Put on the cover (3) and fix it **simultaneously** with the two tension nuts (4).
4. Slip the PVC- water hose over the wide spreading spray diffuser and connect the discharge hose to the outlet of the sieve pan.
5. Fix the holding device and the separate control gear to the wall.

Please notice item 10, Test sieving according to ISO 2591-1!

Connect the mains plug, pay attention to the power supply!

Only connect the control gear delivered with the test sieve shaker to the drive!

10. Test Sieving according to ISO 2591-1

During the sieve analysis, a certain quantity of the test product is divided into fractions by one or several test sieves. The difficulty of the sieve analysis is to determine the separation limit of the separation procedure which is carried out by sieving. It is therefore recommended to determine the analytical separation limit.

Details of the basic principals for sieve analysis are included in the standard sheets.

- **DIN 66 165, part 1 - Basic standard - edition 1987**
- **DIN 66 165, part 2 - Execution - edition 1987**

These sheets are available at the Beuth-Verlag, Berlin.

10.1. Largest particle to be permitted on a test sieve


To avoid damage to the test sieve, the size of the largest particle in the charge should not exceed $10w^{0,7}$ mm, where w is the nominal aperture width in mm.

Example	
Nominal aperture size, w mm	Size of largest particle, apx. mm
0,045	1
0,25	4
1	10
4	25

11. Recommended bulk volume for test sieves with 200 mm \varnothing

1	2	3
Nominal aperture width, w	Bulk volume of material	
μm	Volume of charge app. cm^3	Max. volume of residue app. cm^3
25	22	11
32	26	13
45	30	15
63	35	17
90	42	21
125	50	25
180	60	30
250	70	35
355	80	40
500	100	50
710	120	60
mm		
1	140	70
1,4	160	80
2	200	100
2,8	240	120
4	350	175
5,6	400	200
8	500	250
11,2	800	400
16	1000	500
22,4	1600	800

12. Starting and operating

The sieving time, interval and amplitude can be fed by one pair of entry keys for each function. With the left key the set point **descends**, with the right key it **ascends**. When pressing one of the keys for a longer period, the numbers pass quickly.

Starting the machine

Switch on the machine with the **main switch** (on the bottom side of the Separate Control Gear), all indications of the display read 0.

Zeit / Time

The sieving time is preset in minutes (99 minutes maximum) and appears in the indication field of set values. After the start the preset time runs backwards. The remaining time appears in the indication field of actual values (lower indication field). Once the preset time has run down, the test sieving is finished. Without feeding the sieving time (set point 00) the machine operates continuously.

Interval

Indication 00 = permanent use without interruption (do not use). Indication 01 - 99 = Sieving time in seconds, interruption 1 second. The indication field of actual values shows the interval time running down.

Amplitude

The preset value for the amplitude can be fed by steps of 0.1mm up to 3.0 mm.

The indication field of actual values shows the value of the amplitude measured during the test sieving.

Start / Stop

This key starts the machine. The remaining sieving time appears in the indication field of actual values. During the test sieving process the insert keys for the set points are blocked. After a stop, the remaining sieving time appears in the indication field of actual values and the insert keys for the set points are deblocked.

Continue

By pressing this key a test sieving interrupted with „Stop“ is continued.

Store

By pressing this key 10 set points (time, interval, amplitude) can be stored and **protected from power failure**. Insert the set points and press the key „Store/Speichern“. Insert the storage place (01-10) flashing in the left field of the indication of actual values and press again the key „Store / Speichern“. The values are stored.

Load

By pressing this key the desired storage place (01-10) can be called. After having pressed again on this key and on the key „Start / Stop“, the test sieving can be started.

Delete: Call the storage place by pressing the key „Load/Laden“, set the values to 0 and press the key „Store / Speichern“.

13. Guarantee

All functions of the sieving machine are checked and adjusted before its delivery.

The machine guarantee is for **24 months**.

Test sieves are excluded from guarantee!

We are only liable for defects resulting from flaws or faults in material. The guarantee does not apply to damages caused by transport. Repairs resulting by improper handling or wear due to operation are excluded from guarantee. The guarantee is void, if interventions such as variations or improper repairs are effected which are outside our control.

14. Maintenance and repair

In order to localize defects, the LED-display indicates for a short moment (appx. 5 seconds) a corresponding code of defect as 'E nn' in the first line:

00 - 04

If this indication appears, the machine has an electronic defect and has to be sent back.

Please note that HAVER Test Sieve Shakers only should be repaired, serviced and checked by authorized electro-experts following the national standards (e.f. according to VBG 4). Improper repairs may result in considerable dangers for the user.

If a part of the machine and its equipment is changed or replaced, this part has to be checked again according to EN 60 204-1!

In any case of failure we recommend to contact our service:

HAVER & BOECKER
PARTICLE SIZE ANALYSIS

Ennigerloher Straße 64
59302 OELDE, Germany
Phone: 02522-30363 (GERMANY)
+49-2522-30330 (ABROAD)
Fax: +49-2522-30404
E-Mail: PA@haverboecker.com

The HAVER-Test Sieve Shaker EML 200 digital plus should be checked and serviced at least every 2 years at our works!

15. Disposal of used machine


Used electric and electronic machines contain materials that can be recycled, but also materials that cannot be used again and may cause environmental harm. Do not dispose the used machine with the residual waste, but benefit from our offer to take our used machine back at no charge. It has to be clean and can be delivered carriage free.

16. List of spare parts

EML 200 digital plus T

EML 200 digital plus N

Item No.	Piece	Designation	Article No.
1	2	Nuts with handlebars	00560074
	2	Nuts HAVER TwinNut	00560210
2	1	Inspection glass without hole for machine cover (dry)	00561499
	1	Inspection glass with hole for machine cover (wet)	00561507
3	1	Cover with inspection glass	00561521
	1	Cover without inspection glass	00561514
	1	Cover with inspection glass and wide spreading spray diffuser	00561545
4	3	Rubber buffer	66000074
5		Test sieves	
6	2	Guide rods, clamping system Classic / TwinNut, 445 mm long	00560364
	2	Guide rods, clamping system Classic / TwinNut, 560 mm long	00560009
	2	Guide rods, clamping system Classic / TwinNut, 660 mm long	00560204
7	1	Oscillation magnet	65300016
8	1	Rubber profile for vibration body	00560441
9	1	Vibration body	68500312
10	3	Metal feet rubberized	00560014
11	3	Leaf spring	68500231
12	1	Base plate	68500311
13	1	Separate control unit	00560157
	1	Front panel with control board	00560156
14	1	Rubber ring for vibration body $\varnothing 204 / \varnothing 160$	68500015
15	1	Wide spreading spray diffuser	00561374
	1	Wide spreading spray diffuser with inspection glass	00560489
16	1	Sieve pan with outlet	00540032
	1	Sieve pan without outlet	00540018
17	1	Plug connection with main switch and fuse box	
	1	Plate with cable with screw fitting	
18	2	Counter nut	68500033
19	1	Precision fuse, 20x5 mm, 3,15 Ampere MT (not shown in the drawing)	65110807
20	1	Rubber seal for cover (not shown in the drawing)	00560434

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